

TOBOGA 71-S SPLIT BASKET STRETCHER



SKU: 35200005



PRODUCT DESCRIPTION

The Split Basket Stretcher is made with a high-density polyethylene shell supported by a permanently attached, heavy-duty, aluminum frame and it separates in half for backpacking to a rescue site or for compact storage

Molded-in runners and handholds make it easy to handle, and brass grommets are strategically located for quick use with Model 418-1

Adjustable Lifting Bridle (tag lines should always be used when lifting the basket stretcher). The shell is chemical, UV, rust, and corrosion resistant and a replaceable closed-cell foam pad is attached inside the stretcher for patient comfort.

Four two-piece 5' quick-release patient restraints are simple to use and permit the rescuers to easily change the position of the straps to fit each patient.

Separates in half for backpacking to a rescue site or for compact storage

SPECIFICATIONS

Weight 10 kg

Load limit 272 kg

REXUE LTD. - www.rexue.eu - info@rexue.eu



REXUE
We fire the fire

Width 61 cm
Height 19 cm
Lenght 210 cm



REXUE LTD. - www.rexue.eu - info@rexue.eu